



ORGANIC REACTIONS MECHANISM

335 CHEM - 3

{ Lap Manual }

Prepared by :-

Dr. Amal Fathy

Lec. Sara Hussein

Lec. Eman El-baz

Level 6

1439/1440^H

الرؤية و الرسالة و الأهداف لقسم الكيمياء

Vision, Mission and Goals of Chemistry Department

Vision

To be one of the leading departments in the field of chemistry at the local and global levels.

Mission

The Department of Chemistry, Faculty of Science at Jazan University seeks to develop highly qualified chemists and carry out outstanding scientific research that contributes effectively to the overall development of society by providing a quality educational and research environment that develops knowledge and skills and respects the ethics of education and scientific research.

Goals

- ✚ To provide students with a broad foundation knowledge of the main areas of chemistry: organic, inorganic, analytical, and physical.
- ✚ Students should possess critical thinking and problem-solving abilities.
- ✚ To provide students with the skills required to perform and understand chemical research.
- ✚ To provide students with the skills required to describe, both in writing and orally, chemical processes and procedures.
- ✚ To provide students with the skills required to succeed in a chemical work or related field

الرؤية

أن يكون من الأقسام الرائدة في مجال الكيمياء علي المستوي المحلي والعالمي.

الرسالة

يسعي قسم الكيمياء في كلية العلوم بجامعة جازان لإعداد كيميائيين ذوي كفاءة عالية ، و إنجاز بحوث علمية متميزة تسهم بفاعلية في التنمية الشاملة للمجتمع ، وذلك من خلال توفير بيئة تعليمية و بحثية ذات جودة تنمي المعارف و المهارات و تحترم أخلاقيات التعليم و البحث العلمي.

الأهداف

- ✚ تزويد الطلاب بمعرفة واسعة النطاق من المجالات الرئيسية للكيمياء: العضوية وغير العضوية والتحليلية والفيزيائية والحيوية.
- ✚ يجب أن يمتلك الطلاب قدرات التفكير الناقد وحل المشكلات.
- ✚ تزويد الطلاب بالمهارات اللازمة لأداء وفهم البحوث الكيميائية.
- ✚ تزويد الطلاب بالمهارات المطلوبة لوصف العمليات والإجراءات الكيميائية ، كتابةً أو شفويًا.
- ✚ تزويد الطلاب بالمهارات اللازمة للنجاح في العمل الكيميائي أو مجال آخر ذي صلة.

Safety Rules

- 1- Listen to or read instructions carefully before attempting to do any thing .
- 2- Wear safety goggles to protect your eyes from chemical, hot materials, or things that might be able to shatter .
- 3- Notify your teacher if any spills or accidents occur .
- 4- After handling chemicals , always wash your hands with soap and water .
- 5- During lab work , keep your hands always away from your face .
- 6- Tie back long hair .
- 7- Roll up loose sleeves .
- 8- Know the location of the fire extinguisher , fire blanket , eyewash station, and first aid kit.
- 9- Keep your work area uncluttered . Take to the lab station only what is necessary .
- 10- It is suggested that you wear glasses rather than contact lenses .
- 11- Never put anything into your mouth during a lab experiment .
- 12- Clean up your lab area at the conclusion of the laboratory period .
- 13- Never "horse around" or play practical jokes in the laboratory .

Introduction

Mixture :-

*A mixture contains two or more substances that are chemically not combined.

*Study a mechanism of organic reaction practically involve two phases :

1) Mixture Separation

2) Mixture Components Identification

Organic Compounds :-

1) Acids :- * Aliphatic Acids (Oxalic , Tartaric , Citric Acid).

* Aromatic Acids (Salicylic , Benzoic , Phthalic Acid).

2) Bases :- * Aromatic Amine (*P*-Toluidine , Diphenyl amine , α -Naphthyl amine).

3) Phenols :- * Water soluble phenol (Catechol , Hydroquinone Resorcinol).

*Water insoluble phenol (α - Naphthol , β - Naphthol).

4) Neutral compounds :- * Carbohydrates (Monosaccharaides , Disaccharides ,Starch).

* Salts of Acids (salt of aliphatic acid , salt of aromatic acid).

* Hydrocarbons (Naphthalene , Anthracene).

- The method of separation depends on the solubility of the mixture components in selected solvents.


1 . Mix. {A + A}

Acid + Acid

*Aliphatic Acid + Aromatic Acid

(Notice : All aliphatic acids soluble in cold water, while the aromatic acids insoluble in cold water.)

❖ Separation:-



❖ Identification :- (for aliphatic and aromatic acids)

Preparation of Neutral Solution:- Take a few amount of acid (ppt. , filtrate) in a test tube or beaker. Add NH_4OH gradually drop by drop to the acid till it just alkaline (using litmus paper). Boil the solution until the odor of ammonia is completely removed.

Precipitate :-

1- FeCl_3

N.S of ppt. + FeCl_3

Observation	Result
Violet color	Salicylic Acid
Buff ppt.	Benzoic or phthalic acid

2- Fluorescein Test

(Test to differentiate between benzoic and phthalic acid).

In a clean dry beaker , fuse carefully few crystals of resorcinol and an equal quantity of solid acid (ppt.) then add 2-3drops of $\text{H}_2\text{SO}_4(\text{conc.})$. Pour into excess of NaOH.

**** (If the aromatic acid in form of salt solution add conc. HCl first then filter)**

Observation	Result
Green fluorescence	Phthalic acid
(-ve)No fluorescence	Benzoic acid

Filtrate :- **CaCl_2 Test****N.S filtrate + CaCl_2**

Observation	Result
White ppt. immediately on cold	Oxalic Acid
White ppt. after warming and scratching	Tartaric Acid
White ppt. after boiling	Citric Acid

Conformation Tests

N.S +Amm. $\text{AgNO}_3 \xrightarrow{\Delta}$

Observation	Result
White ppt. $\xrightarrow{\Delta}$ Still white ppt.	Oxalic Acid
white ppt. after $\xrightarrow{\Delta}$ gray $\xrightarrow{\Delta}$ sliver mirror	Tartaric Acid
white ppt. $\xrightarrow{\Delta}$ Gray ppt. only.	Citric Acid


2 . Mix. {A + B}

Acid + Base

*Base + Aliphatic Acid

* Base + Aromatic Acid

❖ Separation:-



❖ **Identification :-** (for aliphatic and aromatic acids it was as illustrated previously at page no.(4,5) ,but it's no need to prepare N.S as the filtrate is in form of salt of acid)


❖ **Chemical properties of Aromatic Amine :-**

1- Azo –dye Test

-Prepare 3 tubes as shown in the figure, allow them to cool in an ice bath, then pour test tube (1) to (2) and finally pour into (3)

(1) $\text{NaNO}_2 + \text{H}_2\text{O}$, (2) solid + dil. HCl , (3) β -naphthol + excess NaOH

Obs.	Res.
(-ve)	Diphenyl Amine
Red orange color	<i>p</i> -toluidine



2-Effect of FeCl_3

- ppt.+dil HCl + FeCl_3

Obs.	Res.
Green color	Diphenyl Amine
(- ve)	<i>p</i> -toluidine

3- (oxidation) Test

- ppt.+dil HNO_3 + drops of conc. H_2SO_4

Obs.	Res.
Dark or blue violet color appear	Diphenyl Amine
(- ve)	<i>p</i> -toluidine


3 . Mix. {A + Ph}

Acid + Phenol



There are four types :-

- (1) Aliphatic Acid + water soluble phenol
- (2) Aliphatic Acid + water insoluble phenol
- (3) Aromatic Acid + water insoluble phenol
- (4) Aromatic Acid + water soluble phenol

❖ Separation:-



- To make sure if the mixture is partially soluble or completely insoluble in water we can do one of two tests.



❖ **Identification :-** (For aliphatic and aromatic acids it was as illustrated previously ,page no.(4,5)

- For type (1) completely soluble (aliphatic acid + water soluble phenol) , directly identify each one in the presence of the other, so no need the step of separation.

Chemical properties of phenol soluble in water :-

1- Effect of FeCl_3

- Place 1 ml of phenol solution in test tube then add few drops of FeCl_3

Obs.	Res.
Green color $\xrightarrow{\text{NaOH}}$ red	Catechol
Red/green crystals then turn to yellow	Hydroquinone
Violet color	Resorcinol

2- Chloroform Test


Phenol solution + chloroform + NaOH (30%) then heat gently(W.B)

Obs.	Res.
Green color	Catechol
Brown color	Hydroquinone

Chemical properties of phenol insoluble in water :-

1- Azo – Dye Test

- Prepare 3 tubes as shown in the figure below , allow them to cool in an ice bath, then pour test tube (1) to (2) and finally pour into (3)



(1) $\text{NaNO}_2 + \text{H}_2\text{O}$, (2) $\text{AnilineHCl} + \text{H}_2\text{O}$, (3) phenol + excess NaOH

Obs.	Res.
Reddish brown ppt.	α - Naphthol
Red ppt.	β -Naphthol

2- Chloroform Test

- In a test tube dissolve phenol in ethanol +chloroform +NaOH then heat gently in(w.B)

Obs.	Res.
Dark blue color turned into green fast.	α - Naphthol & β -Naphthol


4 . Mix. {B +Ph}

Base + Phenol


Two main types :-

- 1) Base + water soluble phenol.
- 2) Base + water insoluble phenol.

❖ Separation:-



- To make sure if the mixture is partially soluble or completely insoluble in water, we can test the filtrate.



Identification :- [for bases and phenols it was as illustrated previously ,page no(7 , 10)].



5 . Mix. {A +N}

Acid + Neutral

(Notice :- Neutral compounds are Carbohydrates , salts of acids and hydrocarbons)

❖ Separation:-

Molisch's Test :- 1 ml of soln. + 1 ml of α -naphthol + conc. H_2SO_4 \longrightarrow violet ring spreads by shaking



Identification :-

* Identify acids as illustrated in page no.(4, 5) . For carbohydrates as follow:-

1-Barfoid Test

1 ml of soln. +2-3ml of Barfoid soln.then heat in W.B

EXP.	Obs.	Res.
Barfoid Test	Red ppt. after 3-5 min. —————>	Monosaccharide
	Red ppt. after 5-7 min. —————>	Disaccharide

Differentiation between Monosaccharides(Glucose – Fructose and Glactose)

2- Rapid furfural

1ml of sugar soln. + 1 ml of α -naphthol + 5-7ml of conc.HCl then heat

EXP.	Obs.	Res.
Rapid furfural	immediate violet color appear —————>	Fructose
	violet color appears after boiling for 1-2min —————>	Glucose
	(-ve) —————>	Glactose

Differentiation between disaccharides (Maltose – Lactose and Sucrose)

3- Fehling's Test


1 ml of sugar soln. + Fehling soln. heat in W.B

EXP.	Obs.	Res.
Fehling's Test	Reddish brown ppt. (Cu_2O) —————>	Maltose – Lactose (which have reduction properties)
	(-ve) —————>	Sucrose (which has no reduction properties)

Differentiation between Maltose and Lactose

4- Osazone Formation

Dissolve 0.2g of sugar +0.4g of phenylhydrazine hydrochloride + 0.6 of sod.acetate in 5 ml water, then heat on W.B for 20min. cool. The osazone separates a yellow crystalline ppt. transfer some of the wet crystals to a microscope slide and try to see the shape of the crystals.

Maltosazone —————> needles shape 

Lactosazone —————> Spherical shape



Mix. {A +N}

Acid + Salt of acid

At first we must determine the type of salt if it's a metallic salt or an ammonium salt

NaOH Test

Mix. + NaOH (30%)

shake well or warm

NH₃ odor


amm.salt of acid

No odor (-ve)

Metallic salt of acid

Types of acid + salt of acid :-

- a) Aliphatic acid + salt of aromatic acid
- b) Aromatic acid + salt of aromatic acid
- c) Aromatic acid + salt of aliphatic acid



[Notice :- salts of acids give all characteristic reactions of the corresponding acids without preparing neutral solution , identify acids as illustrated in page no.(4, 5)].


Mix. {A +N}

Acid + Hydrocarbon


❖ Separation :-

Types of acid + hydrocarbon :-

- a) Aliphatic acid + hydrocarbon
- b) Aromatic acid + hydrocarbon



*In case of the mixture is hydrocarbon+ aliphatic acid, it's better to separate the mix. as follow :-



❖ Identification :-

Chemical properties of hydrocarbons :-

Picrate Test


PPT. + 1-2 drops of picric acid $\xrightarrow{\Delta}$ to dryness then cool

EXP.	Obs.	Res.
Picrate Test	yellow crystals \longrightarrow	Naphthalene
	Red crystals \longrightarrow	Anthracene


6 . Mix. {B +N}

Base + Neutral

Molisch's Test



Base + salt of acid or hydrocarbon




Base +salt of acid

Base +hydrocarbon

Base + carbohydrat


(A) Base + Carbohydrate




*Identify base & sugar as illustrated in page no.(7,13)

(B) Base+ salt of acid

* First the type of salt should be determine (NaOH Test).



(C) Base + hydrocarbon



*Identify base and hydrocarbon as illustrated in page no.(7 , 15)



7 . Mix. {N +N}

Neutral + Neutral


There are four types :-

- (A) Carbohydrate + salt of acid
- (B) Carbohydrate + hydrocarbon
- (C) Salt of aliphatic + salt of Aromatic
- (D) Salt of Acid + hydrocarbon


❖ Separation :-




Type (B) Carbohydrate + hydrocarbon :-



Type (C) Salt of Aliphatic acid + Salt of Aromatic acid : * First the type of salt should be determined (NaOH 30% Test).



Type (D) Salt of acid + Hydrocarbon :-* First the type of salt should be determined (NaOH Test).



*identify sugar & hydrocarbon as illustrated in page no.(13 ,15).